

PORTABLE FREEZER LTFRE34-6



PORTABLE FREEZER LTFRE34-6

LTFRE34-6 PORTABLE FREEZER

Real time temperature value display

Refrigeration and temperature control to ensure the consistency of samples

Easy to operate and could be used repeatedly

Meet the basic requirements of cold chain transportation and storage of vaccines



SPECIFICATIONS

Model	LTFRE34-6
Inner Capacity	33L
Outer size	600x390x320 mm
Product Configuration	1 box 1 digital thermometer 5 95KPa biosafety transport tanks (6 cans are optional) 5 sponge reagent tubes 5 pieces of absorbent cotton 11 crystal ice boxes
Material	PP
Insulation Material	PU
Gross weight	7 kg
Packing size	600x400x330 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com