

PORTABLE FREEZER LTFRE34-7



PORTABLE FREEZER LTFRE34-7

LTFRE34-7 PORTABLE FREEZER

Real time temperature value display

Refrigeration and temperature control to ensure the consistency of samples

Easy to operate and could be used repeatedly

Meet the basic requirements of cold chain transportation and storage of vaccines



SPECIFICATIONS

| | |
|-----------------------|---|
| Model | LTFRE34-7 |
| Inner Capacity | 55L |
| Outer size | 640x395x395 mm |
| Product Configuration | 1 box 1 digital thermometer 8 95KPa biosafety transport tanks (optional 12 tanks) 8 sponge reagent tubes 8 pieces of absorbent cotton 11 crystal ice boxes |
| Material | PP+PE |
| Insulation Material | PU |
| Gross weight | 10 kg |
| Packing size | 650x410x410 mm |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com