

PORTABLE FREEZER LTFRE34-9



PORTABLE FREEZER LTFRE34-9

LTFRE34-9 PORTABLE FREEZER

Real time temperature value display

Refrigeration and temperature control to ensure the consistency of samples

Easy to operate and could be used repeatedly

Meet the basic requirements of cold chain transportation and storage of vaccines



SPECIFICATIONS

Model	LTFRE34-9
Inner Capacity	85L
Outer size	640x520x375 mm
Product Configuration	1 box 1 digital thermometer 20 95KPa biosafety transport tanks 20 sponge reagent tubes 20 pieces of absorbent cotton 15 crystal ice boxes
Material	PP+PE
Insulation Material	PU
Gross weight	15 kg
Packing size	660x540x395 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com