

PORTABLE FREEZER LTFRE35-1



PORTABLE FREEZER LTFRE35-1

LTFRE35-1 PORTABLE FREEZER

AC & DC available Able to used in cars



SPECIFICATIONS

| | |
|--------------------|---|
| Model | LTFRE35-1 |
| Capacity | 30 L |
| Temperature Range | -25 to 30°C |
| Alarm | High&low temp, control error, power failure |
| External Dimension | 610x360x460 mm |
| Power | DC12V; AC100~220V, 50/60Hz |
| Optional Parts | Basket; Datalogger |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com