

PORTABLE FREEZER LTFRE35-2



PORTABLE FREEZER LTFRE35-2

LTFRE35-2 PORTABLE FREEZER



SPECIFICATIONS

Model	LTFRE35-2
Capacity	80 L
Temperature Range	-25 to 30°C
Alarm	High&low temp, control error, power failure
External Dimension	930x580x570 mm
Power	DC12V; AC100~220V, 50/60Hz
Optional Parts	Basket; Datalogger



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com