

STANDARD FUME HOOD LTFUM10-1



STANDARD FUME HOOD LTFUM10-1

Used in Metal analysis, Pharmacy, Research Field, Commercial, Medical, Laboratory.
Also known as Standard Fume Hood, Laboratory Fume Hood, Chemical Fume Hood.

LTFUM10-1 STANDARD FUME HOOD

Microprocessor control system, LED display
Made of porcelain white PP, Resistant to acid, alkali and anti-corrosion
With memory function in case of power-failure
Manual front glass window, height adjustable
Built-in PP centrifugal blower; low noise, easy installation
Chemical resistant phenolic resin work table
Optional: PP work table, epoxy resin board or ceramic board



SPECIFICATIONS

| | |
|-----------------------------|---|
| Model | LTFUM10-1 |
| Air Velocity | 0.3~0.8 m/s |
| Max Opening | 740 mm |
| Work Surface Height | 820 mm |
| Exhaust Duct | 4.3 meters, ϕ 250 mm |
| LED Lamp | 12Wx2 |
| Front Window | Resistant to acid and alkali, Manual, 5 mm toughened glass, height adjustable |
| Blower | Built-in PP centrifugal blowers, speed adjustable |
| Work Table Material | Chemical resistant Phenolic resin |
| Standard Accessory | LED light; Work table: Chemical resistant Phenolic resin; 4.3 meter PVC duct, Diameter: 250 mm pipe strap x2, RVV power cord (3x1.0 mm ² x 3000mm), Base cabinet, Waterproof socket x2, Tail insert, Water tap (PP), Gas tap (PP), Water sink, PP centrifugal blower x 1 |
| Optional Accessory | PP work table, epoxy resin board or ceramic board; Active carbon filter; outer PP centrifugal blower x1, Splicable base cabinet |
| Internal Size (WxDxH) | 787x560x700 mm |
| External Size (WxDxH) | 1047x800x2450 mm |
| Package (Main Body) | 1188x938x1612 mm |
| Package Size (Base Cabinet) | 1188x888x1000 mm or 1200x1000x450 mm |
| Noise Level | \leq 68 dB |
| Gross Weight | 137 kg (Main body) 89 kg or 80 kg (Base Cabinet) |
| Consumption | 400 W |
| Power Supply | AC220V \pm 10%, 60/50Hz; 110V \pm 10%, 60Hz |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com