

FAN FILTER UNIT FFU LTFUM15-1



FAN FILTER UNIT FFU LTFUM15-1

Used in Laboratories, Liquid crystal manufacturing, Research Field, Bio-pharmaceutical.

Also known as Tools for clean air.

LTFUM15-1 FAN FILTER UNIT FFU

FFU shell and HEPA filter adopt split design, easy to install and replace.

Microprocessor control system, energy-saving and low noise.

Long service-life blower, up to 100,000 hours, maintenance-free and small vibration.

Suitable for using separately or connecting together for assembly line.

Made of aluminium coated zinc plate, anti-corrosion, rust-proof and light weight.

SPECIFICATIONS

Model	LTFUM15-1
Air Velocity	0.45m/s±20%
HEPA Filter	99.999% efficiency at 0.3µm
Shell Material	Aluminium coated zinc plate
Blower	Built-in centrifugal blower
Noise	≤60dB
HEPA Filter Size (WxDxH)	970x69x570mm
External Size (WxDxH)	975x575x380mm
Package Size(WxDxH)	1135x540x800mm
Gross Weight	60kg
Consumption	110W
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 50/60Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com