

# STANDARD FUME HOOD LTFUM8-6



# STANDARD FUME HOOD LTFUM8-6

Used in Metal analysis, Pharmacy, Research Field, Commercial, Medical, Laboratory.  
Also known as Standard Fume Hood, Laboratory Fume Hood, Chemical Fume Hood.

## LTFUM8-6 STANDARD FUME HOOD

- Interior body includes melamine board with acid and alkali resistance function
- LCD display
- Exterior body is of cold-roll steel with anti-bacterial powder
- Working chamber is made of solid chemical-resistant material
- UV lamp to sterilize the working area
- Adjustable air velocity with electrical controlled glass
- Adjustable blower speed of a centrifugal fan



## SPECIFICATIONS

Model	LTFUM8-6
Air Velocity	0.3~0.8 m/s
Max Opening	520 mm
Work Surface Height	750 mm
System Exhaust Volume	780m <sup>3</sup> /h
Front Window	Motorized, 5 mm toughened glass; Anti-UV Exterior Body: cold rolled steel with anti-bacteria powder coating.
Work Table Material	Chemical resistant Phenolic resin
LED Lamp	12Wx 1
UV Lamp	20 W Emission of 253.7 nanometers for most efficient decontamination
Blower	Built-in centrifugal blower, Speed adjustable with 9 levels
Standard Accessory	Water tap, Gas tap, Water sink, base stand, UV lampx2, LED lamp, Socket, 4.3 meters exhaust duct, Active carbon filter, foot switch, Waterproof Socket- 2 pcs
Optional Accessory	HEPA filter
Internal Size (WxDxH)	1080x730x745 mm
External Size (WxDxH)	1200x840x2150 mm
Package (Main Body)	1350x1120x1680 mm
Noise Level	≤60 dB
Gross Weight	253 kg
Consumption	400 W
Power Supply	AC220V±10%, 50/60Hz 110V+10%, 60Hz

## ACCESSORIES

Accessory Code	Name
6200812008	Fluorescent lamp

Accessory Code	Name
6200808009	Gas tap

6200812009	UV lamp
6200812010	Base stand
6200812011	Waterproof socket
6200808008	water tap

6200808010	Water sink
6200812012	4 meters PVC exhaust duct
6200808014	Active carbon filter
6200812013	Foot switch



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)