

# AUTOMATIC HYDROGEN GENERATOR LTGAS10-1



# AUTOMATIC HYDROGEN GENERATOR

## LTGAS10-1

Used in Gas quality monitoring, General laboratory, Gas chromatography, TOC, Industries, Environmental, Petrochemical, Pharmaceutical, Cinical, Forensics..

Also known as Laboratory Hydrogen Generator.

## LTGAS10-1 AUTOMATIC HYDROGEN GENERATOR

Automatic control, constant pressure, constant flow

Membrane separation technique, deoxidation and dehydration device, high purity

Easy operation, eliminate the need for expensive gas cylinder

Safe and reliable



### SPECIFICATIONS

Model	LTGAS10-1
Environmental Temperature	0-45°C
Flow Rate	0-1000 ml / min
Output Pressure	0-0.4 Mpa
Pressure Stability	<0.001 Mpa
Dimension	460x370x360mm
Net Weight	20 kg
Consumption Power	400 W
Power Supply	220 V±10%, 50V



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)