

AUTOMATIC HYDROGEN GENERATOR LTGAS10-2



AUTOMATIC HYDROGEN GENERATOR

LTGAS10-2

Used in Gas quality monitoring, General laboratory, Gas chromatography, TOC, Industries, Environmental, Petrochemical, Pharmaceutical, Cinical, Forensics..

Also known as Laboratory Hydrogen Generator.

LTGAS10-2 AUTOMATIC HYDROGEN GENERATOR

Automatic control, constant pressure, constant flow

Membrane separation technique

Compact, frees up valuable space

Easy operation, eliminate the need for expensive gas cylinder

Safe and reliable



SPECIFICATIONS

Model	LTGAS10-2
Hydrogen purity	99.999%
Output Pressure	0-0.3 MPa
Pressure Stability	< 0.001MPa
Dimension	200x140x290 mm
Net Weight	6 kg
Consumption Power	100 W
Power Supply	220V±10%, 50Hz

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com