

AUTOMATIC HYDROGEN GENERATOR LTGAS10-4



AUTOMATIC HYDROGEN GENERATOR

LTGAS10-4

Used in Gas quality monitoring, General laboratory, Gas chromatography, TOC, Industries, Environmental, Petrochemical, Pharmaceutical, Cinical, Forensics..

Also known as Laboratory Hydrogen Generator.

LTGAS10-4 AUTOMATIC HYDROGEN GENERATOR

Automatic control, constant pressure, constant flow

Membrane separation technique

Deoxidation and dehydration device, high purity

Compact, frees up valuable space

Easy operation, eliminate the need for expensive gas cylinder

Safe and reliable



SPECIFICATIONS

Model	LTGAS10-4
Hydrogen purity	99.999%
Flow Rate	0-500 ml / min
Output Pressure	0-0.6 Mpa
Pressure Stability	<0.001 Mpa
Dimension	370x330x180 mm
Net Weight	10 kg
Consumption Power	250 W
Power Supply	220 V±10%, 50 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com