

# AUTOMATIC HYDROGEN GENERATOR LTGAS10-5



# AUTOMATIC HYDROGEN GENERATOR

## LTGAS10-5

Used in Gas quality monitoring, General laboratory, Gas chromatography, TOC, Industries, Environmental, Petrochemical, Pharmaceutical, Cinical, Forensics..

Also known as Laboratory Hydrogen Generator.

## LTGAS10-5 AUTOMATIC HYDROGEN GENERATOR

SPE (Solid Polymer Electrolytes) technique, without lye  
Automatic control, pressure controller, water shortage detector  
Deoxidation and Dehydration device, high purity  
Easy operation and maintenance, long lifetime  
Compact design, Safe and reliable  
CE approved



### SPECIFICATIONS

Model	LTGAS10-5
Hydrogen purity	99.999%
Flow Rate	0-300 ml / min
Output Pressure	0-0.4 Mpa
Pressure Stability	<0.001 Mpa
Dimension	400x360x220 mm
Net Weight	14 kg
Consumption Power	150 W
Power Supply	220 V±10%, 50 Hz



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)