

AUTOMATIC HYDROGEN GENERATOR LTGAS10-6



AUTOMATIC HYDROGEN GENERATOR

LTGAS10-6

Used in Gas quality monitoring, General laboratory, Gas chromatography, TOC, Industries, Environmental, Petrochemical, Pharmaceutical, Cinical, Forensics..

Also known as Laboratory Hydrogen Generator.

LTGAS10-6 AUTOMATIC HYDROGEN GENERATOR

SPE (Solid Polymer Electrolytes) technique, without lye
Automatic control, pressure controller, water shortage detector
Deoxidation and Dehydration device, high purity
Easy operation and maintenance, long lifetime
Compact design, Safe and reliable
CE approved



SPECIFICATIONS

Model	LTGAS10-6
Hydrogen purity	99.999%
Output Pressure	0-0.4 Mpa
Pressure Stability	<0.003 Mpa
Dimension	460x360x350 mm
Net Weight	20 kg
Consumption Power	150 W
Power Supply	220 V±10%, 50 Hz

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com