

AUTOMATIC NITROGEN AIR GENERATOR LTGAS12-2



AUTOMATIC NITROGEN AIR GENERATOR

LTGAS12-2

Used in GC-FID, FPD, NPD, TCD, AED, GC-MS, LC-MS-MS, ICP / NMR, FTIR / IR, ELSD detector, Purging, Ampule Filling, Thermal instruments, Turbo Evaporator (N2 Concentrator), All Analytical Instruments.

Also known as Laboratory Nitrogen Air Generator.

LTGAS12-2 AUTOMATIC NITROGEN AIR GENERATOR

Stainless steel air tank inside

Deoxidization device

Output gas purity is higher

Inside machine use silicon rubber ring (lower content sulfur)

Increased gas output quality

Compact design, easy using, automatically operation, safety



SPECIFICATIONS

Model	LTGAS12-2
Air Pressure	0.4 Mpa
Nitrogen purity	Oxygen content < 3ppm, dew point-56°C
Air Flow Rate	0-2000 ml / min
Dimension	470x260x380 mm
Net Weight	28 kg
Consumption Power	250 W



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com