

# AUTOMATIC ZERO AIR GENERATOR LTGAS8-1



# AUTOMATIC ZERO AIR GENERATOR LTGAS8-1

Used in Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

Also known as Laboratory Air Generator.

## LTGAS8-1 AUTOMATIC ZERO AIR GENERATOR

High precision output pressure

Water automatically drain out

Two levels stainless steel made purifications

Large inside stainless steel air tank

Have two levels pressure stable valve. Output air pressure much stable

Have over pressure protection, low voltage start up, DC, greatly protect inside compressor



## SPECIFICATIONS

Model	LTGAS8-1
Environmental Temperature	0-45°C
Humidity	<85%
Flow Rate	0-5000 ml / min
Output Pressure	0-0.4 Mpa
Pressure Stability	<0.003 Mpa
Noise	<35 dB (A)
Dimension	480x260x350 mm
Net Weight	22 kg
Consumption Power	250 W
Power Supply	220 V±10% V



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)