

AUTOMATIC ZERO AIR GENERATOR LTGAS8-2



AUTOMATIC ZERO AIR GENERATOR LTGAS8-2

Used in Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

Also known as Laboratory Air Generator.

LTGAS8-2 AUTOMATIC ZERO AIR GENERATOR

High precision output pressure

Water automatically drain out

Two levels stainless steel made purifications

Large inside stainless steel air tank

Have two levels pressure stable valve. Output air pressure much stable

Have over pressure protection, low voltage start up, DC, greatly protect inside compressor



SPECIFICATIONS

Model	LTGAS8-2
Output Pressure	0-0.4MPa
Noise	<35dB
Dimension	420x260x350 mm
Net Weight	22 kg
Consumption Power	150 W
Output Flow Rate	0-2000 ml/min
Environmental Temperature	0-45°C
Humidity	≤85%

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com