

AUTOMATIC ZERO AIR GENERATOR LTGAS8-3



AUTOMATIC ZERO AIR GENERATOR LTGAS8-3

Used in Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

Also known as Laboratory Air Generator.

LTGAS8-3 AUTOMATIC ZERO AIR GENERATOR

The hydrocarbon is changed into carbon dioxide and vapor in the catalyzer with Pt-Pd carrier

Increased gas output quality

Compact design, easy using, automatically operation, safety



SPECIFICATIONS

Model	LTGAS8-3
Output air concentration	<0.1PPM
Input air concentration	<100PPM
Maximum Output Pressure	0-0.6 Mpa
Maximum Flow Rate	1000 ml/min
Power Supply	220V±10% 50~60Hz
Consumption Power	150 W
Working Temperature	400-450°C
Environmental Temperature	0-40°C
Humidity	<85%
Dimension	400x360x220 mm
Net Weight	12 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com