

BLUE LIGHT ELECTROPHORESIS MONITOR LTGEL10-1



BLUE LIGHT ELECTROPHORESIS MONITOR LTGEL10-1

Used in Forensic, molecular biology, genetics, microbiology labs, gel electrophoresis instruments.

LTGEL10-1 BLUE LIGHT ELECTROPHORESIS MONITOR



SPECIFICATIONS

Model	LTGEL10-1
Wavelength	470 nm
LED lifetime	≥50000 hours
LED array	Matrix provides single-side lighting
Shiny surface	122 x 35mm
Dimension	206 x57x 20 mm
Weight	0.56kgs
Input power	DC24V 0.75A

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com