

SLIM WHITE ILLUMINATOR LTGEL9-3



SLIM WHITE ILLUMINATOR LTGEL9-3

Used in Forensic, molecular biology, genetics, microbiology labs, gel electrophoresis instruments.

LTGEL9-3 SLIM WHITE ILLUMINATOR

Uniformity-Light guide plate with diffuser design

Easy Operation- Light source is under the gel, avoiding bubble and reflection



SPECIFICATIONS

Model	LTGEL9-3
Light Source	White light LEDs
Viewing Size	240Dx160W mm
Dimensions (DxWxH)	292x217x12.8 mm
Power	12V DC
Weight	1.2 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com