

# PATHOLOGY WORKSTATION LTHEP10-2



# PATHOLOGY WORKSTATION LTHEP10-2

Used in Hospital pathology department, Pathology laboratory.

Also known as Pathological sampling bench.

## LTHEP10-2 PATHOLOGY WORKSTATION

Odor-remove function

Single water sink, two-person operations

Exhaust air system: High-power suction device, low noise

Optional water heater is available, realizing hot & cold water adjustment

The base cabinet is inward-concave type, convenient for sitting operation

Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes

Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning



## SPECIFICATIONS

Model	LTHEP10-2
Work Surface Height	850 mm
Exhaust Airflow Volume	500~600m <sup>3</sup> /h
Main Body Material	304 stainless steel
Work Table Material	304 stainless steel
Package Size (WxDxH)	1660x890x2170 mm
Gross Weight	300 kg
Working Area Size (WxDxH)	1395x430x900 mm
External Size (WxDxH)	1500x750x1950 mm
Consumption	300 W
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)