

UV STERILIZATION CABINET LTHEP13-1



UV STERILIZATION CABINET LTHEP13-1

Used in Hospital , Laboratory , House.
Also known as UV Disinfection Cabinet.

LTHEP13-1 UV STERILIZATION CABINET



SPECIFICATIONS

Model	LTHEP13-1
Type	With Ozone
UV Ray Wavelength	253.7/184.9nm
Sterilization Way	2 sides of UV ray+1 side of Ozone
Capacity	200 L
Power	120 W
UV Lamp Material	Quartz glass
Shelves	Standard 2pc, can be optional
Internal Size	518x523x778 mm
External Size	750x663x1200 mm
Packing Size	875x788x1265 mm
Gross Weight	100 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com