

UV STERILIZATION CABINET LTHEP13-2



UV STERILIZATION CABINET LTHEP13-2

Used in Hospital , Laboratory , House.
Also known as UV Disinfection Cabinet.

LTHEP13-2 UV STERILIZATION CABINET



SPECIFICATIONS

Model	LTHEP13-2
Type	Without Ozone
UV Ray Wavelength	253.7nm
Sterilization Way	3 sides of UV ray
Capacity	200 L
Power	120 W
UV Lamp Material	Quartz glass
Shelves	Standard 2pc, can be optional
Internal Size	518x523x778 mm
External Size	750x663x1200 mm
Packing Size	875x788x1265 mm
Gross Weight	100 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com