

AIR PURIFIER LTHEP14-1



AIR PURIFIER LTHEP14-1

Used in Hospital, small clinic, laboratory, office, meeting room and houses..

Also known as Air Purifier.

LTHEP14-1 AIR PURIFIER

Low noise

LCD display

Remote control

Adjustable air velocity

Thorough sterilizing with UV Lamp

Mobile Type, Wall Hanging Type, Suction Ceiling Type



SPECIFICATIONS

Model	LTHEP14-1
Filter	Primary Filter, HEPA filter
HEPA Filter	99.99% efficiency at 0.3pm
UV Lamp	Germicidal UV lamp, Emission of 253.7 nanometers for most efficient decontamination
Optimal Space	≤20m ²
Noise	≤58 dB
Consumption	240 W
Adjustable Air Speed	3 levels
Cleanness At Outlet	Class 100
Gross Weight	78 kg
Package Size (WxDxH)	1800x440x1050 mm
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz
External Size (WxDxH)	900x120x1500 mm

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com