

TISSUE PROCESSOR LTHIS10-2



TISSUE PROCESSOR LTHIS10-2

Used in Clinical and Research Histopathology.

Also known as Tissue Embedding Cassettes.

LTHIS10-2 TISSUE PROCESSOR

It has a capacity of 200 cassettes, set a new standard for specimen quality and reliability in automated tissue processing. designed for routine clinical and research histopathology, is an innovative, smart processor for paraffin infiltration of tissue. Straight forward routine user operations and a variety of 'Smart' features improves specimen quality and laboratory economy.



SPECIFICATIONS

Model	LTHIS10-2
Capacity	200 cassettes 165kg
Nominal voltage	AC220V
Nominal frequency	60 Hz
Main fuses	10 A
Nominal power	1000 VA
Dry weight, unpacked	approx. 100 kg
Weight, packed	130 kgs/ 165kgs/172kgs
Operating temperature range	15 °C to 35 °C
Relative humidity	10 % to 80 %, non-condensing



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com