

DOUBLE BASKET TISSUE PROCESSOR LTHIS10-6



DOUBLE BASKET TISSUE PROCESSOR LTHIS10-6

Used in Clinical and Research Histopathology.

Also known as Tissue Embedding Cassettes.

LTHIS10-6 DOUBLE BASKET TISSUE PROCESSOR

This machine can run both single hangingbasket and double basket, one dehydrator is equivalent to two hanging baskets. Large LCD touch screen, all English menu, simple and clear operation. Adopt intelligent program control technology to protect equipment and specimens. Intelligent separate protection specimen in case of failure. When the dehydrator fails before the anhydrous alcohol cylinder, the tissue is protected in the fixed cylinder. when the dehydrator fails before the paraffin cylinder and after the anhydrous alcohol cylinder, the tissue is protected in the anhydrous alcohol cylinder; when the hanging basket fails in the paraffin cylinder dehydrator, the tissue is protected in the paraffin cylinder. Adopt new water heating mode, low pressure contactless temperature control technology, more safe. Reagent cylinder can be fully heated to avoid the effect of dehydration due to the difference of temperature in winter and summer.



SPECIFICATIONS

Model	LTHIS10-6
Volume of individual reagent containers	2000 ml
Number of reagent containers	14
Hanging basket filtrate	Time can be adjusted arbitrarily
Number of dehydrated tissues	More than 200
Processing time in a single container	0 ~ 99 hours adjustable
Power supply	AC 220 V \pm 10 % 50 Hz
Power	\leq 1500 W



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com