

AUTOMATIC TISSUE PROCESSOR LTHIS10-8



AUTOMATIC TISSUE PROCESSOR LTHIS10-8

Used in Clinical and Research Histopathology.

Also known as Tissue Embedding Cassettes.

LTHIS10-8 AUTOMATIC TISSUE PROCESSOR

Based on different tissues, several running programs available for customers. Double-control heat without water. Self-provided power supply can ensure to continue 6 hours' working even if the electricity has been cut off. With transmission of gear and rack, mechanical failures can be reduced. Several standby programs available



Product Image Coming Soon

SPECIFICATIONS

Model	LTHIS10-8
Tanks	12 Tanks (3 wax Tanks in the end)
Capacity of single tank	1500 ml
Basket can hold	75 specimens
Single-tanks processing time	Randomly set within the range from 0 to 99 hour(s).
The temperature of wax cylinder	Can be randomly set within the range of 99°C
Time for cylinder to be lifted	38 seconds
Tissue-protection cylinder sequence	7th cylinder
Times of mixing per minute	3 times



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com