

AUTOMATIC TISSUE PROCESSOR LTHIS10-9



AUTOMATIC TISSUE PROCESSOR LTHIS10-9

Used in Clinical and Research Histopathology.

Also known as Tissue Embedding Cassettes.

LTHIS10-9 AUTOMATIC TISSUE PROCESSOR

Prompt judge and automatic processing in case of abnormal condition, which can prevent the tissue from being damaged.

Double-control heat without water. Based on different tissues, several running programs available for customers. With transmission of gear and rack, mechanical failures can be reduced. Several standby programs available.



SPECIFICATIONS

Model	LTHIS10-9
Tanks	12 Tanks (3 wax cylinders in the end)
Capacity of single- tanks	1000 ml
Basket can hold	60 specimens
Single tanks processing time	set arbitrarily in the range of 0 to 99 hour(s)



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com