

IMMUNO AUTOSTAINER LTHIS11-5



IMMUNO AUTOSTAINER LTHIS11-5

Used in Histology, Cytology, Hematology, Pathology, and Microbiology staining.

Also known as Tissue dyeing, Slide Stainer.

LTHIS11-5 IMMUNO AUTOSTAINER

Open system, Suitable for various kinds of specimen. True hands free, walk away ,overnight operation.Easy Windows™ based programming,password protection. Allows for several stainers using single computer simultaneously. User-specified protocols/User's own reagents,several runs per day. Robust, proven construction and operation. No disposable pipette tips. Per-run reagent volume checks.



SPECIFICATIONS

Model	LTHIS11-5
Slides Capacity	48
Reagents Capacity	64
Minimum dropping reagent	50 µl
Dropping reagent increment	50 µl, 100 µl, 150 µl, 200 µl
Maximal output power	≤ 350 VA
Working noise	≤ 60 dB (A)
Application	IHC,HE more
Dimension (mm)	1015x665x520
Power supply	220 V or 110 V
Rated current	1.6 A
Weight	116 kgs



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com