

CASSETTE PRINTER LTHIS12-2



CASSETTE PRINTER LTHIS12-2

Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens.
Also known as Cassettes printer.

LTHIS12-2 CASSETTE PRINTER

Laser printing and easy for maintaining. Capability of linking with the management system of the pathology department. Print online or off-line. The writing is clear and delicate, never wears out or fades, which meets the characteristic requirements of keeping slides for long-term reference. No need for preheating, can be immediately used after starting up, safe and noiseless, and free of Maintenance.



SPECIFICATIONS

Model	LTHIS12-2
Capacity	150 Slides
Average printing speed	≤ 3S / slide
Printing accuracy	≥ 1250 pi
Maximum Capacity	≥ 150 slides
Minimum Capacity	≥ 7 slides
Ambient Humidity	45 – 85 %
Dimensions	465x200x361 mm (LengthxWidthxHeight)
Laser Power	Pulse laser ≤ 2 W
System interface	Integrate with LIS/HIS/P IS
Printing Methods	Print online or off-line
Loading method	Cartridge clip type
External device	Support Bar-code scanner inputting
Consumable items	5×75 mm, 1 - 1.2 mm glass slide



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com