

PATHOLOGICAL SLIDE SCANNER LTHIS12-6



PATHOLOGICAL SLIDE SCANNER LTHIS12-6

Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens.

Also known as Cassettes printer.

LTHIS12-6 PATHOLOGICAL SLIDE SCANNER

The PrimeHisto XE is a histological slide scanner. With easy manipulation/setting, PrimeHisto XE allows the histological samples to be digitized in a matter of seconds. Ability to be suitable for most histological slide samples and can scan the entire histological slide image in a scanning area of 24.3x36.5mm. Also, has a high resolution of up to 10,000 DPI.



Product Image Coming Soon

SPECIFICATIONS

Model	LTHIS12-6
Resolution	10000x10000 DPI
Dynamic range	3.8 OD
Scan area	24.3x36.5 mm
Scanning medium	75x25 mm section
Photo formats	TIFF, JPEG
Light	transmits LED white light
Operating systems	Windows XP/Vista / 7/8/10, Mac OS X 10.5+
Volume	260x155x70 mm
Weight	1.05 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com