

PARAFFIN DISPENSER LTHIS13-4



PARAFFIN DISPENSER LTHIS13-4

Used in Pathology, Histology.

Also known as Tissue Embedding Station.

LTHIS13-4 PARAFFIN DISPENSER

The full automatic procedure control & new heating element gives feature of quick heating & energy saving. It possesses the function of memorization to keep the set temperature automatically after startup. The LED shows heating state.



SPECIFICATIONS

Model	LTHIS13-4
Range of temperature	Ambient~100 °C
Precision	± 1 %
Voltage	AC 110 V 60 Hz
Power	1200 VA
Capacity	10000 ml
Overall dimension	350x450x510 mm
Net Weight	17 kg

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com