

PARAFFIN TRIMMER LTHIS13-5



PARAFFIN TRIMMER LTHIS13-5

Used in Pathology, Histology.

Also known as Tissue Embedding Station.

LTHIS13-5 PARAFFIN TRIMMER

The grooved heating surface channels melted wax to a removable collection tray. The collection tray connects magnetically to the base of the main unit. Disposable plastic liners are included to insert into the collection tray for easy disposal of runoff wax. This unit now includes a digital control panel for personalized settings.



SPECIFICATIONS

Model	LTHIS13-5
Working area dimension	150x220 mm
Temperature range	Ambient to 75 °C
Temperature control accuracy	± 1 °C
Supply voltage	220 V ± 10 % 50 / 60 Hz OR 110 V ± 10 % 50 / 60 Hz
Max power	240 W
Outside dimensions (WxDxH)	220x160x175 mm
Weight	2.5 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com