

TISSUE WATER BATH

LTHIS13-6



TISSUE WATER BATH LTHIS13-6

Used in Pathology, Histology.
Also known as Tissue Embedding Station.

LTHIS13-6 TISSUE WATER BATH

Heat conduct fast, durable. Self memory of temperature. Display presetting temperature & current temperature. Corrosion-resistance of the bath. With the function of spreading and baking slices, the cylinder can be taken out, and 9 glass slides can be placed on the baking slices.



SPECIFICATIONS

Model	LTHIS13-6
Range of setting temperature	Ambient to 90 °C (constant adjusting)
Temperature precise error	+ 2 °C
Allowance ambient temperature	0 – 40 °C
Working voltage	220 V / 50 Hz
Power	250 W
Dimensions overall	450x400x200 mm (WxDxH)
Weight	7 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com