

TISSUE BAKING PROCESSOR LTHIS13-7



TISSUE BAKING PROCESSOR LTHIS13-7

Used in Pathology, Histology.

Also known as Tissue Embedding Station.

LTHIS13-7 TISSUE BAKING PROCESSOR

Heat conduct fast, durable. Self memory of temperature. Display presetting temperature & current temperature. Corrosion-resistance of the bath.

SPECIFICATIONS

Model	LTHIS13-7
Range of setting temperature	Ambient to 90 °C (constant adjusting).
Temperature precise error	+ 2 °C.
Allowance ambient temperature	0 – 40 °C
Working voltage	220 V / 50 Hz
Power	250 W
Dimensions overall	450x400x200 mm (WxDxH)
Weight	7 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com