

PORTABLE ULTRASONIC PROCESSOR LTHMG11-1



PORTABLE ULTRASONIC PROCESSOR

LTHMG11-1

Used in Biotechnology, Pharmaceutical, Chemical, Industrial, Environmental.

Also known as Portable Ultrasonic Homogenizer Sonicator, Portable ultrasonic cell crusher.

LTHMG11-1 PORTABLE ULTRASONIC PROCESSOR



SPECIFICATIONS

Model	LTHMG11-1
Operating frequency range	20-25 KHz, frequency automatic tracking
The end of the horn diameter	Φ 3 or 5
Ultrasonic Power	150W, can be adjusted, hit into ways to fine-tune each of 5W
Parameters sample handling capacity	0.1-100ml (need to select the appropriate probe)
Generator output power	150W (max)
Operating voltage	220V±5%



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com