

# PORTABLE ULTRASONIC PROCESSOR LTHMG11-2



# PORTABLE ULTRASONIC PROCESSOR

## LTHMG11-2

Used in Biotechnology, Pharmaceutical, Chemical, Industrial, Environmental.

Also known as Portable Ultrasonic Homogenizer Sonicator, Portable ultrasonic cell crusher.

## LTHMG11-2 PORTABLE ULTRASONIC PROCESSOR



### SPECIFICATIONS

Model	LTHMG11-2
Operating frequency range	20-25 KHz, frequency automatic tracking
The end of the horn diameter	Φ 3 or 5
Ultrasonic Power	250W, can be adjusted, hit into ways to fine-tune each of 5W
Parameters sample handling capacity	0.1-100ml (need to select the appropriate probe)
Generator output power	250W (max)
Operating voltage	220V±5%



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)