

MICROPLATE SEALER LTHMG8-1



MICROPLATE SEALER LTHMG8-1

Used in Pharmaceutical, medical, Food Industry, Beverage Industry.

Also known as Capping Machine.

LTHMG8-1 MICROPLATE SEALER

Easy to switch the cap opening module, one equipment can be achieved 96 holes rack or 100 holes rack cap opening and closing operation.

Quick to switch the cap hat and suitable for all kinds of brand sample tube cap opening and closing.

Built-in integration with intelligence sensor which can be identified different volume sample tube.

It not only can be independent operation but also can be together with other liquid workstation and robotic arm comprise a automatic system.

Drawer type inlets and outlets for avoiding sample contamination.

Touch screen control panel One-button operation provide the user with direct and convenient operation feeling.



SPECIFICATIONS

Model	LTHMG8-1
Rack type available	96 or 100 Holes
Screen	Touch Screen control panel
Sensor	Built-in integration with intelligence sensor



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com