

PLATE SEALER LTHMG8-2



PLATE SEALER LTHMG8-2

Used in Pharmaceuticals, Mobile inspection laboratories, Disease prevention and control centers, Food Industry, Beverage Industry.

Also known as Screw on capping machine.

LTHMG8-2 PLATE SEALER

One-button start

Simple and fast operation

Safe and labor-saving

Capping/sealing speed $\leq 4s$



SPECIFICATIONS

Model	LTHMG8-2
Throughput	8-96, suitable for 8-strip PCR tubes and 96-well PCR plates
Start Method	Button
Speed	Capping\sealing speed $\leq 4s$
Function	8-strip PCR tubes capping, 96-well PCR plate capping, PCR plate sealing film
Adapter	2 types, 8-strip PCR tubes and 96-well PCR plates, can be customized
Pressure	Peak pressure: 300kg
External Size	277x194x230 mm
Package Size	300x260x340 mm
Weight	Net Weight- 9.2kg Gross Weight- 12kg
Power Supply	100-240V 50/60Hz

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com