

MICROPLATE SEALER LTHMG8-3



MICROPLATE SEALER LTHMG8-3

Used in Pharmaceuticals, Mobile inspection laboratories, Disease prevention and control centers, Food Industry, Beverage Industry.

Also known as Screw on capping machine.

LTHMG8-3 MICROPLATE SEALER

Fully automatic and one key operation, no need of additional manual operation

Adapts to varieties of plate types

Adjustable sealing pressure

A variety of communication interfaces are reserved, which can be quickly connected to other machines

Each roll of film can seal more than 350 plates(reference quantity)



SPECIFICATIONS

Model	LTHMG8-3
Positions	1pc
Method	Cold Seal
Speed	<16s/plate
Operating Mode	Manual plate change, automatic sealing
Application	PCR plate, deep well plate and micro-plate conforming to SBS standard(Customizable)
Control Mode	Buttons
Interface	USB
External Size	310x380x410 mm
Package Size	454x382x464 mm
Weight	Net Weight- 16.5kg Gross Weight- 20.5kg
75 W	Power Supply



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com