

BARCODE READER LTHMG9-1



BARCODE READER LTHMG9-1

Used in Laboratory, Research, Medical.

Also known as Tube Reader, Tube Scanner, Tube Barcode scanner, Sample Scanner.

LTHMG9-1 BARCODE READER

It can be finish within 5 seconds for one of 96 holes rack or 384 holes rack.


Advanced decoder technology can be identified stable and secure even if the code is fuzzy.

Optional item for 1-D code rack scanning function.

It can be scanned difference of height sample tube in the same rack .

Compact, ergonomic design assures efficient of workflow.

Detection of rack type automatically for setting optimum scanning.



SPECIFICATIONS

Model	LTHMG9-1
Rack type available	13,24,48,96,384 holes rack
1-D Code available	2/5 Industrial / Interleaved, Code 39, Code 128, Pharmacode, Codabar, EAN 13
2-D Code available	Datamatrix ECC 200, PDF417, QR Code
Camera	10 mega-pixel camera CMOS



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com