

ICE FLAKER LTICE10-10



ICE FLAKER LTICE10-10

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

LTICE10-10 ICE FLAKER

- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin



SPECIFICATIONS

Model	LTICE10-10
Capacity	200 kg/24 h
Storage Bin Capacity	60 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	680x611x1100 mm
Packing Dimension	750x675x1150 mm
Weight	80/85 kg
Power	1040 W
Power Supply	220 V, 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com