

ICE FLAKER LTICE10-11



ICE FLAKER LTICE10-11

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

LTICE10-11 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



SPECIFICATIONS

Model	LTICE10-11
Capacity	300 kg/24 h
Storage Bin Capacity	150 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	560x800x1700 mm
Packing Dimension	650x850x1750 mm
Weight	95/105 kg
Power	1360 W
Power Supply	220 V, 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com