

BULLET ICE MAKER LTICE10-12



BULLET ICE MAKER LTICE10-12

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

LTICE10-12 BULLET ICE MAKER

- Transparent Bullet Ice
- Air Cooling
- Stainless Steel 304 Exterior
- Adjustable Feet
- Slide Away Access Door
- World Famous Brand Compressor
- Fully Automatic Controller
- Improved Design For Higher Reliability



SPECIFICATIONS

Model	LTICE10-12
Capacity	120 kg/24 h
Storage Bin Capacity	55 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	680x611x1050 mm
Packing Dimension	750x675x1150 mm
Weight	70/75 kg
Power	800 W
Power Supply	220 V /110 V, 50 / 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com