

CUBE ICE MAKER LTICE10-15



CUBE ICE MAKER LTICE10-15

Used in Healthcare, Food industries, Refrigerated perishables, Biochemical and biological studies, Medical laboratories, Beverage industries, Pharmaceutical.

Also known as Ice Cuber Machine, Laboratory Cube Ice Maker.

LTICE10-15 CUBE ICE MAKER

High grade stainless steel SUS304 is used which is anticorrosive and durable. The compact structure is simple and space-saving.

Fully automatically computer controlled, convenient for use.

Automatic detection and protection in case of operation, water deficiency, ice full and malfunction.

Highly efficient fluoride-free compressor, energy saving and environment-friendly.

Suitable for Pubs, KTV and dancing halls, gyms, cold drinks shops, milk tea shops, cafes, convenience stores, restaurants and other entertainment venues.

Split-type structure, large ice bin, large ice storage capacity and strong applicability.



SPECIFICATIONS

Model	LTICE10-15
Capacity	500 kg/24 h
Storage Bin Capacity	280 kg
Refrigerant	R134a / R404a
Condenser Type	Water cooled
Machine Head	760x600x740 mm
Cabinet	760x800x980 mm
Packing Machine Head	845x705x950 mm
Packing Cabinet	800x850x1060 mm
Weight	123/133 kg
Power Supply	220 V, 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com