

# CUBE ICE MAKER LTICE10-16



# CUBE ICE MAKER LTICE10-16

Used in Healthcare, Food industries, Refrigerated perishables, Biochemical and biological studies, Medical laboratories, Beverage industries, Pharmaceutical.

Also known as Ice Cuber Machine, Laboratory Cube Ice Maker.

## LTICE10-16 CUBE ICE MAKER

High grade stainless steel SUS304 is used which is anticorrosive and durable. The compact structure is simple and space-saving.

Fully automatically computer controlled, convenient for use.

Automatic detection and protection in case of operation, water deficiency, ice full and malfunction.

Highly efficient fluoride-free compressor, energy saving and environment-friendly.

Suitable for Pubs, KTV and dancing halls, gyms, cold drinks shops, milk tea shops, cafes, convenience stores, restaurants and other entertainment venues.

Split-type structure, large ice bin, large ice storage capacity and strong applicability.



## SPECIFICATIONS

Model	LTICE10-16
Capacity	1000 kg/24 h
Storage Bin Capacity	205 kg
Refrigerant	R134a / R404a
Machine Head	1230x630x740 mm
Cabinet	1230x920x1095 mm
Packing Machine Head	1320x680x950 mm
Packing Cabinet	1290x1000x1160 mm
Weight	205/225 kg
Power Supply	380 V, 60 Hz



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)