

# ICE FLAKER LTICE10-2



# ICE FLAKER LTICE10-2

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

## LTICE10-2 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



## SPECIFICATIONS

|                      |                |
|----------------------|----------------|
| Model                | LTICE10-2      |
| Capacity             | 30 kg/24 h     |
| Storage Bin Capacity | 10 kg          |
| Refrigerant          | R134a          |
| Condenser Type       | Air cooled     |
| Overall Dimension    | 330x470x605 mm |
| Packing Dimension    | 405x555x665 mm |
| Weight               | 28/31 kg       |
| Power                | 300 W          |
| Power Supply         | 220 V, 60 Hz   |

## OPTIONAL ACCESSORIES

| Accessory Code | Name         |
|----------------|--------------|
| 2900808008     | Water filter |



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)