

# ICE FLAKER LTICE10-2



# ICE FLAKER LTICE10-2

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

## LTICE10-2 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



## SPECIFICATIONS

Model	LTICE10-2
Capacity	30 kg/24 h
Storage Bin Capacity	10 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	330x470x605 mm
Packing Dimension	405x555x665 mm
Weight	28/31 kg
Power	300 W
Power Supply	220 V, 60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name
2900808008	Water filter



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)