

ICE FLAKER LTICE10-4



ICE FLAKER LTICE10-4

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

LTICE10-4 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



SPECIFICATIONS

| | |
|----------------------|----------------|
| Model | LTICE10-4 |
| Capacity | 50 kg/24 h |
| Storage Bin Capacity | 15 kg |
| Refrigerant | R134a |
| Condenser Type | Air cooled |
| Overall Dimension | 400x543x720 mm |
| Packing Dimension | 470x615x780 mm |
| Weight | 36/39 kg |
| Power | 380 W |
| Power Supply | 220 V, 60 Hz |

OPTIONAL ACCESSORIES

| Accessory Code | Name |
|----------------|--------------|
| 2900808008 | Water filter |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com