

ICE FLAKER LTICE10-7



ICE FLAKER LTICE10-7

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

LTICE10-7 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



SPECIFICATIONS

Model	LTICE10-7
Capacity	85 kg/24 h
Storage Bin Capacity	25 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	400x510x825 mm
Packing Dimension	470x590x875 mm
Weight	45/49 kg
Power	480 W
Power Supply	220 V, 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2900808008	Water filter



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com