

ICE FLAKER LTICE10-8



ICE FLAKER LTICE10-8

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

LTICE10-8 ICE FLAKER

- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin



SPECIFICATIONS

Model	LTICE10-8
Capacity	100 kg/24 h
Storage Bin Capacity	40 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	500x611x940 mm
Packing Dimension	570x685x1000 mm
Weight	58/63 kg
Power	520 W
Power Supply	220 V, 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2900808008	Water filter



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com